


ELECTROPHORESIS MANUAL SYSTEM

PICOPHORI



Seleo S.r.l.
I Traversa Bugnano sc
81030 Orta di Atella (CE) Italy
Tel. +39 (0)81 5022877 - (0)81 2359553
P.IVA 08850061212
Homepage: www.seleo.it
e-mail: info@seleo.it / sales1@seleo.it

MADE IN
ITALY

ELECTROPHORESIS MANUAL SYSTEM

The system consists of:

Electrophoretic Power Supply

Electrophoretic Chamber

Multiple Work Support with Application Plate



Samples plate at 40 positions for Serum Proteins and Hemoglobins

Applicator for Serum proteins and Hemoglobins

Applicator for Immunofixation

Samples plate at 5 positions for Immunofixation

PICOPHOR



Electrophoretic Power Supply

PicoPhor electrophoretic power supply unit can be programmed both in Voltage and in Current thanks to its multifunction selector, its display and its easy-to-use menu system.

Immunofixation

picoPhor offers all the simplicity and speed possible in preparing samples of proteins Immunofixation.

Best paired with microPhor densitometer

Complete your system with our microPhor densitometer for reading the results of Serum Proteins and Hemoglobins electrophoresis



SELEO SRL

Your partner in Protein Electrophoresis