

 **YDRA 09**

Highly Reliable Photometer



SELEO SRL

MADE IN
ITALY



Highly Reliable Photometer

SELEO reagents preprogrammed.

Smart card technology ensures reliable results with SELEO consumables and validated reagents.

Guided workflow on touchscreen LCD

FEATURES

Hydra 09 is a photometer that allows the execution of turbidimetric and chromogenic methods.

The possibility of using pre-calibrated kits with memory cards allows the running of non-routine exams.

By using a touchscreen interface and an intuitive manager program, Hydra 09 is very simple to use and reliable as routine or back up photometer. Available preprogrammed.

EXAMPLES OF TESTS

- | | | | |
|----------------|---------------------|-----------------|--------------|
| ● Homocysteine | ● Protein C | ● Valproic Acid | ● FOB |
| ● HbA1C | ● Protein S | ● Carbamazepine | ● Calprotein |
| ● G-6-PDH | ● Kappa free chain | ● Phenobarbital | ● D-Dimer |
| ● Ammonium | ● Lambda free chain | ● Digoxina | ● Fibrogen |

TECHNICAL FEATURES

- | | |
|---|---|
| ● Halogen Lamp | ● Q.C. features |
| ● Up to 10 optical filters | ● Visual control during work phases |
| ● 8 sample position in pre-incubation | ● USB port for data transfer or software upgrade |
| ● 1 reading cell | ● RJ45 (ethernet) port for interface LIS |
| ● Minimum reaction volume of 200 ul | ● Languages: Italian, English (more languages on request) |
| ● Memory card with calibration curves | ● LCD color touch screen 7" easy to use with an innovative user friendly software |
| ● Open or pre-programmed system | ● Dimensions: 34x22x10cm |
| ● More than 100.000 results stored | ● Weight: 3 Kg |
| ● Integrated thermal printer | ● Remote assistance via AnyDesk |
| ● Temperature controlled by electrical thermostat at 37° +/-0,1 | |
| ● Methods: Kinetic; End-point; Fixed-time | |



Prima Traversa Bugnano s.c.
81030 Orta di Atella (CE) Italy
Tel. +39 081 5022877-2359553
sales1@seleo.it - info@seleo.it
www.seleo.it